



78529

10X Phosphate Buffered Saline (PBS) for
molecular biology

Part B

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
pH (10X Buffer)	6.9 ± 0.1
pH (1X Buffer)	7.4 ± 0.1
Sterility (Bacterial count)	Nil
DNase, RNase, Protease	Not detected

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	36 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 ml	38229090 (GST 12%)
500 ml	38229090 (GST 12%)

Available Packages

100 ml
500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com