



78586

**L-Phenylalanine Methyl Ester Hydrochloride
extrapure, 99%**

Part B

CAS : 7524-50-7

Molecular Formula : C10H13NO2.HCl

Molecular Weight : 215.69

Specifications

Appearance (Colour)	White to off - white
Appearance (Form)	Crystalline powder
Assay (ex Cl)	min. 99%
Melting Point	158 - 160°C
Specific rotation [?]	+37° to +38° (C=2, Ethanol)

General Information

Storage	0 to 4°C (Freeze)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	29379020 (GST 18%)
25 Gms	29379020 (GST 18%)
5 Gms	29379020 (GST 18%)
Type of Packing	
1 Gms	Glass Bottle
5 Gms	Glass Bottle

Available Packages

1 Gms
5 Gms
25 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com