



79007

## Nitromethane pure, 97%

Part A

CAS : 75-52-5

Molecular Formula : CH<sub>3</sub>NO<sub>2</sub>

Molecular Weight : 61.04

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to pale yellow, oily
Appearance (Form)	Liquid
Colour (APHA)	max. 30
Solubility with water	Slightly soluble
Assay (GC)	min. 97%
Density (g/ml) @ 20°C	1.132-1.134
Refractive Index (20°C)	1.381-1.383
Boiling Range	100-102°C
Water (KF)	max. 0.5%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
HSN Code	
500 ml	29042090 (GST 18%)
2500 ml	29042090 (GST 18%)
25 Ltr	29042090 (GST 18%)

### Available Packages

500 ml

2500 ml

25 Ltr

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com