



79167

Paraffin Liquid Heavy (Mineral Oil Heavy) extrapure, 110-230mPas

Part A

CAS : 8012-95-1, 8042-47-5

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Viscous liquid
Density (g/ml) @ 20°C	0.860 - 0.904 g/ml
Viscosity (20°C)	110 - 230 mPas
Solid paraffin	Passes
Acidity	Passes test
Sulphur Compound (S)	Passes
Readily carbonisable substances	Passes

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
25 Ltr	27129090 (GST 18%)
2500 ml	27129090 (GST 18%)
500 ml	27129090 (GST 18%)

Available Packages

500 ml

2500 ml

25 Ltr

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com