



79563

Phenolphthalein Indicator 1% Soln in Ethanol

Part A

CAS : 77-09-8

Specifications

| | |
|--------------------------|----------------------------------|
| Appearance (Clarity) | Clear |
| Appearance (Colour) | Colourless to slight pale yellow |
| Appearance (Form) | Liquid |
| Suitability as indicator | Passes test |

General Information

| | |
|---------------------|---------------------------------------|
| Storage | 25 to 40°C (Room Temperature) |
| Shelf Life | 24 Months |
| IMDG Identification | Not Regulated for Transport (Non-Haz) |
| HSN Code | |
| 125 ml | 38229090 (GST 12%) |
| 500 ml | 38229090 (GST 12%) |

Available Packages

125 ml
500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com