



79576

Olive Oil

Part A

CAS : 8001-25-0

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Pale yellow to yellow to greenish yellow
Appearance (Form)	Liquid
Refractive Index (20°C)	1.468 - 1.470
Acid value	max. 0.3
Water (KF)	max. 2%
Suitability as enzyme Substrate	Passes test

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 ml	15099090 (GST 5%)
250 ml	15099090 (GST 5%)

Available Packages

100 ml

250 ml

500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com