



79862

## Stannic Chloride 1M in Dichloromethane

Part A

CAS : 7646-78-8

Molecular Formula :  $\text{SnCl}_4$

Molecular Weight: 260.52

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Pale yellow to yellow
Appearance (Form)	Liquid
Concentration	0.95 - 1.10 M

### General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	60 Months
IMDG Identification	UN No.:1827 , IMCO Class No.:8 , Packing Group:II
HSN Code	
100 ml	38229090 (GST 12%)
500 ml	38229090 (GST 12%)

### Available Packages

100 ml

---

500 ml

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com