



79897

3-(Hydroxymethyl) Phenylboronic acid extrapure, 95%

Part B

CAS : 87199-15-3

Molecular Formula : C₇H₉BO₃

Molecular Weight : 152.00

Specifications

| | |
|---------------------|--------------------|
| Appearance (Colour) | White to off-white |
| Appearance (Form) | Crystalline powder |
| Assay (NT) | min. 95.0% |

General Information

| | |
|------------|-------------------------------|
| Storage | 25 to 40°C (Room Temperature) |
| Shelf Life | 60 Months |
| HSN Code | |
| 1 Gms | 29319090 (GST 18%) |
| 10 Gms | 29319090 (GST 18%) |
| 5 Gms | 29319090 (GST 18%) |

Available Packages

1 Gms

5 Gms

10 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com