



80053

Bismuth Metal Lumps, 99.5%

Part A

CAS : 7440-69-9

Molecular Formula : Bi

Atomic Weight : 208.98

Specifications

|                        |            |
|------------------------|------------|
| Appearance (Colour)    | Shiny grey |
| Appearance (Form)      | Lumps      |
| Assay (EDTA Titration) | min. 99.5% |

General Information

|                     |                                       |
|---------------------|---------------------------------------|
| Storage             | 25 to 40°C (Room Temperature)         |
| Shelf Life          | 60 Months                             |
| IMDG Identification | Not Regulated for Transport (Non-Haz) |
| HSN Code            |                                       |
| 100 Gms             | 81061020 (GST 18%)                    |
| 500 Gms             | 81061020 (GST 18%)                    |
| Type of Packing     |                                       |
| 100 Gms             | Plastic Bottle                        |
| 500 Gms             | Plastic Bottle                        |

Available Packages

|         |
|---------|
| 100 Gms |
| 500 Gms |

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com