



80145

DAPI Dihydrochloride (4,6-Diamidino-2-Phenylindole Dihydrochloride) Solution (1mg/ml in H₂O)

Part B

CAS : 28718-90-3

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	slight pale yellow
Appearance (Form)	Liquid

Other Information

Application	For nucleic acid staining.
-------------	----------------------------

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	12 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 ml	38229090 (GST 12%)
5 x1 ml	38229090 (GST 12%)
Type of Packing	
1 ml	Glass Bottle

Available Packages

1 ml
5 x1 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com