



80233

Potassium Dichromate extrapure AR, 99.9%

Part A

CAS : 7778-50-9

Molecular Formula : K₂Cr₂O₇

Molecular Weight : 294.18

Specifications

Appearance (Colour)	Orange
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Orange
Assay (Iodometric)	min. 99.9%
Loss on drying	max. 0.05%
Chloride (Cl)	max. 0.001%
Sulphate (SO ₄)	max. 0.005%
Iron (Fe)	max. 0.002%
Heavy Metals (Pb)	max. 0.005%
Copper (Cu)	max. 0.001%
Sodium (Na)	max. 0.02%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3086 , IMCO Class No.:6.1(5.1) , Packing Group:I
HSN Code	
5 Kg	28415090 (GST 18%)
500 Gms	28415090 (GST 18%)
Type of Packing	
500 Gms	Plastic Bottle

Available Packages

500 Gms

5 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com