



80400

## Papain 5x USP ex. Papaya Latex, 30000USP U/mg

Part B

CAS : 9001-73-4

### Specifications

Appearance (Colour)	Off-white to Pale brown
Appearance (Form)	Powder
Activity	min. 30000 USP units/mg
Unit Definition	One USP unit of Papain activity is that which releases the equivalent 1 µg of tyrosine from a specified casein substrate under the conditions of the assay, using the enzyme concentration that liberates 40 µg of tyrosine per mL of sample solution.
Enzyme solution for activity determination	Prepared in 0.06 M phosphate buffer, pH 6.0.

### General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	35079099 (GST 18%)
25 Gms	35079099 (GST 18%)
500 Gms	35079099 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle

### Available Packages

25 Gms  
100 Gms  
500 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com