



80726

## Sodium Alginate LV extrapure, 4-12cps, 30mesh

Part A

CAS : 9005-38-3

### Specifications

Appearance (Colour)	Pale yellowish brown
Appearance (Form)	Powder
pH (1% aq. solution)	6.0 - 8.0
Viscosity (1% aq. at 25°C)	4 - 12cps
Loss on drying	max. 15%
Gelation Test	Passes test
Sulphated Ash	30 - 36%
Mesh size	30mesh

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
5 Kg	39131010 (GST 18%)
500 Gms	39131010 (GST 18%)
Type of Packing	
500 Gms	Plastic Bottle

### Available Packages

500 Gms

5 Kg

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com