



80909

## Tetrapropylammonium Chloride extrapure, 98%

Part A

CAS : 5810-42-4

Molecular Formula : C<sub>12</sub>H<sub>28</sub>ClN

Molecular Weight : 221.81

### Specifications

Appearance (Colour)	White to off - white
Appearance (Form)	Crystalline compound
Solubility (Turbidity) 10% aq. solution	Clear
Assay (On dry basis)	min. 98%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
10 Gms	29239000 (GST 18%)
50 Gms	29239000 (GST 18%)
Type of Packing	
50 Gms	Glass Bottle

### Available Packages

10 Gms

50 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com