



81085

Triethyl Orthoformate extrapure, 99%

Part A

CAS : 122-51-0

Molecular Formula : C₇H₁₆O₃

Molecular Weight : 148.20

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Assay (GC)	min. 99%
Density (g/ml) @ 20°C	0.890 - 0.892
Refractive Index (20°C)	1.391-1.393
Boiling Range	143-146°C
Peroxide (H ₂ O ₂)	max. 0.05%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2524 , IMCO Class No.:3 , Packing Group:III
HSN Code	
2500 ml	29151300 (GST 18%)
500 ml	29151300 (GST 18%)
Type of Packing	
2500 ml	Glass Bottle
500 ml	Glass Bottle

Available Packages

500 ml

2500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com