



81499

Sodium Tungstate Dihydrate extrapure, 98%

Part A

CAS : 10213-10-2

Molecular Formula : Na₂WO₄.2H₂O

Molecular Weight : 329.86

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% hot aq. solution	Clear
Solubility (Colour) 10% hot aq. solution	Colourless
Assay (T)	min. 98%
pH (2% aq. solution)	7.0 - 10.0
Chloride (Cl)	max. 0.05%
Sulphate (SO ₄)	max. 0.05%
Nitrate (NO ₃)	max. 0.01%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	28418090 (GST 18%)
2.5 Kg	28418090 (GST 18%)
250 Gms	28418090 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle

Available Packages

100 Gms

250 Gms

2.5 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com