



81618

Silicon Antifoaming Agent 30% (Non-ionic)

Part A

Specifications

Appearance (Colour)	White
Appearance (Form)	Thick fine suspension
pH	6.5-7.5 (as is)

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 ml	39100090 (GST 18%)
500 ml	39100090 (GST 18%)
Type of Packing	
100 ml	Plastic Bottle
500 ml	Plastic Bottle

Available Packages

100 ml

500 ml

1000 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com