



81658

Acetanilide extrapure, 99%

Part A

CAS : 103-84-4

Molecular Formula : C8H9NO

Molecular Weight : 135.17

Specifications	
Appearance (Colour)	White to off-white
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 5% solution in Methanol	Clear to slightly opalescent
Assay	min. 99%
Melting Point	114 - 116°C
Sulphated Ash	max. 0.05%

General Information	
Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
500 Gms	29242910 (GST 18%)
Type of Packing	
500 Gms	Plastic Bottle

Available Packages	
500 Gms	

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com