



81862

Seralose® CL-2B (2% Beaded Agarose Cross-Linked Susp., 60-200micron)

Part B

Specifications

Appearance (Form)	Pre-swollen 2% Beaded agarose aqueous suspension crosslinked with Epichlorohydrin
Flow rate	min. 100 mL/hr (5.0 cm seralose bead in a 1.5 cm diameter column with 20 cm water head.)
Particle Size	60 - 200 microns
Matrix	2% Agarose
Stability pH range	3.0 - 14.0
Fractionation range for proteins	7 x 10 ⁴ - 4 x 10 ⁷
Fractionation range for polysaccharides	1 x 10 ⁵ - 2 x 10 ⁷
DNA Exclusion limits	1352 base pairs
Approx. Exclusion limits for proteins	45 x 10 ⁶
Approx. Exclusion limits for polysaccharides	20 x 10 ⁶
Chromatography technique	Affinity coupling
Chemical stability	Stable to commonly aqueous buffers
Autoclavable	Autoclavable 20 min. at 120°C in pH 7

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
HSN Code	
100 ml	13023100 (GST 18%)
500 ml	13023100 (GST 18%)

Available Packages

100 ml
500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com