



82045

## Silicone Oil 330-370 CPS

Part A

CAS : 63148-62-9

Molecular Formula : [-Si(CH<sub>3</sub>)<sub>2</sub>O-]<sub>n</sub>

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Almost colourless
Appearance (Form)	Liquid
Refractive Index (20°C)	1.401 - 1.405
Density (g/ml) @ 20°C	0.966 - 0.977
Viscosity (20°C)	330 - 370 cP
Boiling Point	>250°C

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1000 ml	39100020 (GST 18%)
250 ml	39100020 (GST 18%)
2500 ml	39100020 (GST 18%)

### Available Packages

250 ml

---

1000 ml

---

2500 ml

---

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com