



82149

Benzaldehyde extrapure AR, ExiPlus, Multi-Compendial, 99%

Part A

CAS : 100-52-7

Molecular Formula : C₇H₆O

Molecular Weight : 106.12

Specifications

ExiPlus	meets compendial specs of BP,Ph Eur
Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Assay (GC)	min. 99%
Density (g/ml) @ 20°C	1.043 - 1.045
Refractive Index (20°C)	1.544 - 1.545
Boiling Range	177 - 179°C
Acidity (C ₆ H ₅ COOH)	max. 0.3%
Sulphated Ash	max. 0.01%
Ionised Chloride (Cl)	max. 0.0001%
Water (KF)	max. 0.1%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1990 , IMCO Class No.:9 , Packing Group:III
HSN Code	
2500 ml	29122100 (GST 18%)
500 ml	29122100 (GST 18%)
Type of Packing	
500 ml	Glass Bottle

Available Packages

500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com