



82395

Trimethyltetradecylammonium Chloride extrapure, 99%

Part A

CAS : 4574-04-3

Molecular Formula : $\text{CH}_3(\text{CH}_2)_{13}\text{N}(\text{Cl})(\text{CH}_3)_3$

Molecular Weight : 291.94

Specifications

Appearance (Colour)	White to off white
Appearance (Form)	Crystalline Powder
Assay (NT) (On dry basis)	min.99%

Other Information

Application	Suitable for Ion Pairing & HPLC applications
-------------	--

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29239000 (GST 18%)
25 Gms	29239000 (GST 18%)

Available Packages

25 Gms

100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com