



82523

Celite 545, size 20-40 micrometer

Part A

CAS : 91053-39-3

Specifications

Appearance (Colour)	White to off - white
Appearance (Form)	Powder
pH (10% aq. suspension)	min. 8.5
Particle Size	20 - 40 μ m
Loss on ignition	max. 0.3%

Other Information

Notes	Used as a filter aid.
-------	-----------------------

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Kg	25120090 (GST 5%)
500 Gms	25120090 (GST 5%)
Type of Packing	
1 Kg	Plastic Jar
500 Gms	Plastic Bottle

Available Packages

500 Gms

1 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com