



82634

## MES6: ICP Non-Metals

Part B

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to Pale yellow
Appearance (Form)	Liquid
Arsenic (As)	100.0 ± 0.5 µg/mL
Boron (B)	100.0 ± 0.5 µg/mL
Phosphorus (P)	100.0 ± 0.5 µg/mL
Sulphur (S)	100.0 ± 0.5 µg/mL
Selenium (Se)?	100.0 ± 0.5 µg/mL
Silicon (Si)	100.0 ± 0.5 µg/mL
Tellurium (Te)	100.0 ± 0.5 µg/mL

### General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	36 Months
HSN Code	
100 ml	38229090 (GST 12%)

### Available Packages

100 ml

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com