



82736

1,3-Dioxolane extrapure AR, 99%

Part A

CAS : 646-06-0

Molecular Formula : C₃H₆O₂

Molecular Weight : 74.06

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Assay (GC)	min. 99%
Density (g/ml) @ 20°C	1.060 - 1.070
Boiling Range	74 - 76°C
Refractive Index (20°C)	1.395 - 1.405

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1166 , IMCO Class No.:3 , Packing Group:II
HSN Code	
2500 ml	29329990 (GST 18%)
500 ml	29329990 (GST 18%)

Available Packages

500 ml

2500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com