



82745

## Acetonitrile (ACN) extrapure AR, ACS, ExiPlus, Multi-Compendial, 99.5%

Part A

CAS : 75-05-8

Molecular Formula : CH<sub>3</sub>CN

Molecular Weight : 41.05

### Specifications

ExiPlus	meets compendial specs of USP, BP, Ph.Eur
Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Colour (APHA)	max. 10
Assay (GC)	min. 99.5%
Density (g/ml) @ 20°C	0.780 - 0.783
Refractive Index (20°C)	1.343 - 1.344
Boiling Range	81 - 83°C
Non Volatile Matter	max. 0.001%
Acidity (CH <sub>3</sub> COOH)	max. 0.01%
Iron (Fe)	max. 0.00001%
Heavy Metals (Pb)	max. 0.00001%
Alkalinity	max. 0.001%(NH <sub>3</sub> )
Benzene (GC)	max. 0.1%
Cyanide (CN)	max. 0.005%
Miscibility with water	Passes test
Water (KF)	max. 0.1%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1648 , IMCO Class No.:3 , Packing Group:II
HSN Code	
2500 ml	29269000 (GST 18%)
500 ml	29269000 (GST 18%)
Type of Packing	
2500 ml	Glass Bottle
500 ml	Glass Bottle

## Available Packages

500 ml

2500 ml

---

### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: [info@srchem.com](mailto:info@srchem.com), website [www.srchem.com](http://www.srchem.com)

---