



82818

## Trifluoroacetic Acid (TFA) for HPLC & UV Spectroscopy, 99.9%

Part A

CAS : 76-05-1

Molecular Formula : C<sub>2</sub>HF<sub>3</sub>O<sub>2</sub>

Molecular Weight : 114.02

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless, fuming
Appearance (Form)	Liquid
Assay (T)	min. 99.9%
Density (g/ml) @ 20°C	1.492-1.494
Water (KF)	max. 0.05%
Absorbance of sample in 1cm cell v/s H <sub>2</sub> O	
@260nm	max. 0.9
@270nm	max. 0.1
@280nm	max. 0.05
@290nm	max. 0.04
@300nm	max. 0.03
@320nm	max. 0.025

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2699 , IMCO Class No.:8 , Packing Group:I
HSN Code	
100 ml	29159093 (GST 18%)
25 ml	29159093 (GST 18%)
500 ml	29159093 (GST 18%)
Type of Packing	
25 ml	Glass Bottle
500 ml	Glass Bottle

## Available Packages

25 ml

100 ml

500 ml

---

### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: [info@srchem.com](mailto:info@srchem.com), website [www.srchem.com](http://www.srchem.com)

---