



82863

Tetrahydrofuran (THF) (Stabilized) pure, 99.5%

Part A

CAS : 109-99-9

Molecular Formula : C₄H₈O

Molecular Weight : 72.11

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Assay (GC)	min 99.5%
Density (g/ml) @ 20°C	0.886-0.888
Refractive Index (20°C)	1.406-1.408
Boiling Range	65-67°C
Non Volatile Matter	max. 0.005%
Peroxide (H ₂ O ₂)	max. 0.01%
Acidity (CH ₃ COOH)	max. 0.002%
Water (KF)	max. 0.1%
Stabilizer	0.025% di-tert butyl-p-cresol

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2056 , IMCO Class No.:3 , Packing Group:II
HSN Code	
25 Ltr	29321100 (GST 18%)
2500 ml	29321100 (GST 18%)
500 ml	29321100 (GST 18%)
Type of Packing	
2500 ml	Glass Bottle

Available Packages

500 ml
2500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com