



82949

Benedict's Reagent Quantitative

Part A

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Blue
Appearance (Form)	Liquid
Suitability as reagent for quantitative analysis of Reducing sugar	Passes test

Other Information

Notes

Copper sulphate 18gms ,Sodium carbonate 200gms ,Sodium citrate 200gms ,Potassium thiocyanate 125gms , Potassium ferrocyanide 0.25gms dilute in 1litre water

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	12 Months
IMDG Identification	UN No.:3082 , IMCO Class No.:9 , Packing Group:III
HSN Code	
2500 ml	38229090 (GST 12%)
500 ml	38229090 (GST 12%)

Available Packages

500 ml

2500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com