



83147

## Aluminium Ammonium Sulphate Dodecahydrate extrapure AR, 99.5%

Part A

CAS : 7784-26-1

Molecular Formula :  $\text{AlNH}_4(\text{SO}_4)_2 \cdot 12\text{H}_2\text{O}$

Molecular Weight : 453.32

### Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless
Assay (Complexometric)	min. 99.5%
pH (5% aq. solution)	~ 3.0
Chloride (Cl)	max. 0.001%
Iron (Fe)	max. 0.0005%
Heavy Metals (Pb)	max. 0.0005%
Phosphate (PO <sub>4</sub> )	max. 0.001%
Copper (Cu)	max. 0.0005%
Silicate (SiO <sub>2</sub> )	max. 0.001%
Potassium (K)	max. 0.02%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
500 Gms	28333010 (GST 18%)
Type of Packing	
500 Gms	Plastic Bottle

### Available Packages

500 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com