



83407

## Sodium Citrate Tribasic Dihydrate extrapure, 98%

Part A

CAS : 6132-04-3

Molecular Formula : C<sub>6</sub>H<sub>5</sub>Na<sub>3</sub>O<sub>7</sub>·2H<sub>2</sub>O

Molecular Weight : 294.10

### Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless
Assay (NT)	min. 98%
Acidity	max. 1.0 ml N%
Chloride (Cl)	max. 0.01%
Sulphate (SO <sub>4</sub> )	max. 0.02%
Heavy Metals (Pb)	max. 0.001%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Kg	29181520 (GST 18%)
5 Kg	29181520 (GST 18%)
500 Gms	29181520 (GST 18%)
Type of Packing	
1 Kg	Plastic Bottle
500 Gms	Plastic Bottle

### Available Packages

500 Gms

1 Kg

5 Kg

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com