



83417

Sodium Phosphate Dibasic Dodecahydrate extrapure AR, 99%

Part A

CAS : 10039-32-4

Molecular Formula : $\text{Na}_2\text{HPO}_4 \cdot 12\text{H}_2\text{O}$,

Molecular Weight : 358.14

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless
Assay (HCl titration)	min. 99%
pH (5% aq. solution)	9.0 - 9.6
Chloride (Cl)	max. 0.005%
Sulphate (SO ₄)	max. 0.05%
Iron (Fe)	max. 0.005%
Heavy Metals (Pb)	max. 0.002%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
5 Kg	28352200 (GST 18%)
500 Gms	28352200 (GST 18%)

Available Packages

500 Gms

5 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com