



83438

Charcoal Activated (280) extrapure AR, ACS, ExiPlus, Multi-Compendial

Part A

CAS : 7440-44-0

Molecular Formula : C

Molecular Weight : 12.01

Specifications

ExiPlus	meets compendial specs of USP, BP, Ph.Eur
Appearance (Colour)	Black
Appearance (Form)	Fine powder
pH (2 % aq. suspension)	6.0 - 9.0
MB value	min. 280 mg/g
Loss on drying	max. 10%
Sulphated Ash	max. 4%
Chloride (Cl)	max. 0.1%
Iron (Fe)	max. 0.05%
Sulphate (SO4)	max. 0.05%

Other Information

Notes (Forbidden by Air)

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1362 , IMCO Class No.:4.2 , Packing Group:III
HSN Code	
100 Gms	38029019 (GST 18%)
5 Kg	38029019 (GST 18%)
500 Gms	38029019 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle
5 Kg	Plastic Drum
500 Gms	Plastic Bottle

Available Packages

100 Gms

500 Gms

5 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com
