



83740

Gelatin Powder ex. Porcine (Type B), 120g Bloom

Part A

CAS : 9000-70-8

Specifications

Appearance (Colour)	Pale yellow to yellowish brown
Appearance (Form)	Granular powder
Solubility (Turbidity) 5% hot aq. solution	Clear to slightly hazy
pH (5% aq. solution)	4.0 - 6.0
Loss on drying	max. 15%
Gel Bloom test	120 to 150 gm

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
5 Kg	35030090 (GST 18%)
500 Gms	35030090 (GST 18%)
Type of Packing	
500 Gms	Plastic Bottle

Available Packages

500 Gms

5 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com