



83834

N,N-Dimethyl-L-Phenylalanine extrapure, 98%

Part B

CAS : 17469-89-5

Molecular Formula : C₁₁H₁₅O₂N

Molecular Weight : 193.24

Specifications

Appearance (Colour)	White to Off-white
Appearance (Form)	Powder
Assay (NT)	min. 98.0%
Water (KF)	max. 1.0%
Identification (IR Spectrum)	Passes test

General Information

Storage	25 to 40°C (Room Temperature)
HSN Code	
1 Gms	29224990 (GST 18%)
5 Gms	29224990 (GST 18%)
10 Gms	29224990 (GST 18%)

Available Packages

1 Gms
5 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com