



84198

N,N,N,N-Tetramethyl Ethylenediamine (TEMED) extrapure AR, ExiPlus, Multi-Compendial, 99%

Part A

CAS : 110-18-9

Molecular Formula : C₆H₁₆N₂

Molecular Weight : 116.21

Specifications

ExiPlus	meets compendial specs of Ph.Eur
Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Assay (GC)	min. 99%
Density (g/ml) @ 20°C	0.773 - 0.777 g
Refractive Index (20°C)	1.417 - 1.419
Boiling Range	120 - 122°C
Water (KF)	max. 1%

Other Information

Notes	Electrophoresis and immunoelectrophoresis reagent.
-------	--

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2372 , IMCO Class No.:3 , Packing Group:II
HSN Code	
100 ml	29212990 (GST 18%)
500 ml	29212990 (GST 18%)

Available Packages

100 ml
500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com