



84210

## Silica Gel Fumed 200 (Hydrophilic)

Part A

CAS : 112945-52-5

### Specifications

Appearance (Colour)	White
Appearance (Form)	Very light fine powder
pH (4% aq. suspension)	3.7 - 5.5
Loss on drying	max. 7%
Specific surface area (m <sup>2</sup> /g)	175 - 225 m <sup>2</sup> /g

### Other Information

Notes	Equivalent to Aerosil 200
-------	---------------------------

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	28399090 (GST 18%)
2.5 Kg	28399090 (GST 18%)
250 Gms	28399090 (GST 18%)

### Available Packages

100 Gms

---

250 Gms

---

2.5 Kg

---

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com