



84444

Tetrabutylammonium Dihydrogen Phosphate,  
0.1M in HPLC Water

Part A

CAS : 5574-97-0

Molecular Formula : $\text{C}_{16}\text{H}_{38}\text{NO}_4\text{P}$

Molecular Weight: 339.45

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Molarity (M)	0.095-0.105
Absorbance of 0.005M in water, 1cm cell vs HPLC grade water	
@200nm	max. 0.050
@210nm	max. 0.030
@254nm	max. 0.020
@300nm	max. 0.010
@500nm	max. 0.010

General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 ml	38229090 (GST 12%)
500 ml	38229090 (GST 12%)
Type of Packing	
100 ml	Glass Bottle

Available Packages

100 ml
500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com