



84627

Cadmium Acetate Dihydrate extrapure AR, 99%

Part A

CAS : 5743-04-4

Molecular Formula : $Cd(CH_3COO)_2 \cdot 2H_2O$

Molecular Weight : 266.53

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline compound
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless
Assay	min. 99%
pH (5% aq. solution)	6.7 - 7.3
Chloride (Cl)	max. 0.001%
Sulphate (SO ₄)	max. 0.005%
Iron (Fe)	max. 0.0005%
Heavy Metals (Pb)	max. 0.0005%
Copper (Cu)	max. 0.0005%
Sodium (Na)	max. 0.002%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3077 , IMCO Class No.:9 , Packing Group:III
HSN Code	
100 Gms	29152990 (GST 18%)
500 Gms	29152990 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle
500 Gms	Plastic Bottle

Available Packages

100 Gms

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com