



84697

Indium Metal Ingots, 99.9%

Part B

CAS : 7440-74-6

Atomic Formula : In

Atomic Weight : 114.82

Specifications

Appearance (Form)	Shining metal ingots
Assay	min. 99.9%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Restricted for export from India
HSN Code	
10 Gms	81129900 (GST 18%)

Available Packages

10 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com