



84753

## Methyl Eugenol extrapure, 99%

Part B

CAS : 93-15-2

Molecular Formula : C<sub>11</sub>H<sub>14</sub>O<sub>2</sub>

Molecular Weight : 178.23

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to pale yellow
Appearance (Form)	Liquid
Assay (GC)	min. 99%
Density (g/ml) @ 20°C	1.035 - 1.036 g/ml
Refractive Index (20°C)	1.534 - 1.535

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 ml	29095090 (GST 18%)
2500 ml	29095090 (GST 18%)
500 ml	29095090 (GST 18%)

### Available Packages

100 ml

500 ml

2500 ml

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com