



84905

General Purpose Kit - II (Pre-Weighed Media)

Part D

Other Information

Applications

Nutrient agar cultivates a wide variety of microorganisms, MacConkey agar aids in isolation of gram negative enteric bacteria and Blood Agar Base is used in detecting hemolytic of fastidious microorganisms.

Composition

(63971) Nutrient Agar	28.00 gm/lt.
(76875) MacConkey Agar	50.00 gm/lt.
(84905) Blood Agar Base	40.00 gm/lt.

General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	36 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Kit	38229090 (GST 12%)

Available Packages

1 Kit

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com