



84910

Adenosine-2',3' Cyclic Monophosphate Sodium Salt extrapure, 98%

Part B

CAS : 37063-35-7

Molecular Formula : C₁₀H₁₁N₅O₆PNa

Molecular Weight : 351.10

Specifications

| | |
|--------------------------------|--------------------|
| Appearance (Colour) | White |
| Appearance (Form) | Crystalline powder |
| Assay | min. 98% |
| Spectral Ratio (A250nm/A260nm) | 0.76 - 0.82 |
| Spectral Ratio (A280nm/A260nm) | 0.13 - 0.17 |

General Information

| | |
|---------------------|---------------------------------------|
| Storage | -20 °C (Blue/Dry Ice) |
| Shelf Life | 60 Months |
| IMDG Identification | Not Regulated for Transport (Non-Haz) |
| HSN Code | |
| 100 Mg | 29349990 (GST 18%) |
| 25 Mg | 29349990 (GST 18%) |

Available Packages

25 Mg

100 Mg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com