



85091

Cholesterol Oxidase ex. Streptomyces Sp., 15U/mg solids

Part B

CAS : 9028-76-6

Specifications

Appearance (Colour)	Pale Yellow
Appearance (Form)	Crystalline powder
Activity	min. 15 units/mg powder
Unit Definition	One unit will convert 1.0 micromole of cholesterol to 4-cholesten-3-one per minute at pH 7.5 at 25 °C. Prepared in 20 mM potassium phosphate buffer containing 0.2% BSA, pH 7.0.
Enzyme solution for activity determination	

General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
2500 Units	35079099 (GST 18%)
500 Units	35079099 (GST 18%)

Available Packages

500 Units
2500 Units

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com