



85108

Ethyl Methanesulphonate (EMS) extrapure, 99%

Part B

CAS : 62-50-0

Molecular Formula : C₃H₈SO₃

Molecular Weight : 124.15

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Assay (GC)	min. 99%
Density (g/ml) @ 20°C	1.205 - 1.207 g/ml
Refractive Index (20°C)	1.418 - 1.420

General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	60 Months
IMDG Identification	UN No.:2810 , IMCO Class No.:6.1 , Packing Group:II
HSN Code	
10 Gms	98020000 (GST 18%)
100 Gms	98020000 (GST 18%)

Available Packages

10 Gms
100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com