



85124

Sodium Tripolyphosphate Anhydrous (STPP) extrapure, 57-59% as P2O5

Part A

CAS : 7758-29-4

Molecular Formula : Na₅P₃O₁₀

Molecular Weight : 367.86

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline compound
Assay (P ₂ O ₅ basis)	57 - 59%
pH (1% aq. solution)	9.1 - 10.1
Loss on drying	max. 1%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1000 Gms	28353100 (GST 18%)
5000 Gms	28353100 (GST 18%)
Type of Packing	
1000 Gms	Plastic Bottle
5000 Gms	Plastic Jar

Available Packages

1 Kg

5 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com