



85471

Xylene Dried, 99.5%, 0.005% water

Part A

CAS : 1330-20-7

Molecular Formula : C₈H₁₀

Molecular Weight : 106.17

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Total Assay (GC, C ₈ H ₁₀ isomers + Ethylbenzene)	min. 99.5%
Density (g/ml) @ 20°C	0.860-0.865
Refractive Index (20°C)	1.497-1.499
Boiling Range	136-140°C
Non Volatile Matter	max. 0.001%
Acidity (HCl)	max. 0.001%
Substance Discoloured by H ₂ SO ₄	Passes test
Water	max. 0.005%

Other Information

Notes

Packed with Dr. Dried Seal which are designed with corrosion free Crown Cap & with specially developed coloured septa. This enable repetitive needle puncturing easy and securing the consistency of the product from outside atmosphere.

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1307 , IMCO Class No.:3 , Packing Group:III
HSN Code	
1000 ml	29024400 (GST 18%)
Type of Packing	
1000 ml	Glass Bottle

Available Packages

1000 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com
