



85833

Graphite Nanopowder (Type 1), 400nm, 98%

Part C

CAS : 7782-42-5

Molecular Formula : C

Molecular Weight : 12.01

Specifications

| | |
|---------------------|----------|
| Appearance (Colour) | Black |
| Appearance (Form) | Powder |
| Assay | min. 98% |
| APS | 400nm |

General Information

| | |
|---------------------|---------------------------------------|
| Storage | 25 to 40°C (Room Temperature) |
| Shelf Life | 60 Months |
| IMDG Identification | Not Regulated for Transport (Non-Haz) |
| HSN Code | |
| 100 Gms | 25049090 (GST 5%) |
| 25 Gms | 25049090 (GST 5%) |
| 5 Gms | 25049090 (GST 5%) |
| Type of Packing | |
| 100 Gms | Plastic Bottle |

Available Packages

5 Gms
25 Gms
100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com