



85994

## n-Propanol Dried, 99.5%, 0.005% water

Part A

CAS : 71-23-8

Molecular Formula : C<sub>3</sub>H<sub>8</sub>O

Molecular Weight : 60.10

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Assay (GC)	min. 99.5%
Density (g/ml) @ 20°C	0.803-0.805
Refractive Index (20°C)	1.385-1.386
Boiling Range	96-98°C
Non Volatile Matter	max. 0.002%
Acidity (C <sub>3</sub> H <sub>6</sub> O <sub>2</sub> )	max. 0.001%
Water (KF)	max. 0.005%

### Other Information

Notes

Packed with Dr. Dried Seal which are designed with corrosion free Crown Cap & with specially developed coloured septa. This enable repetitive needle puncturing easy and securing the consistency of the product from outside atmosphere.

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1274 , IMCO Class No.:3 , Packing Group:II
HSN Code	
1000 ml	29051210 (GST 18%)

### Available Packages

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com