



86017

2-Heptanone (Methyl-n-Amyl Ketone, MAK) extrapure 99%

Part B

CAS : 110-43-0

Molecular Formula : $\text{CH}_3(\text{CH}_2)_4\text{COCH}_3$

Molecular Weight : 114.19

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Assay (GC)	min. 99%
Density (g/ml) @ 20°C	0.813-0.816
Refractive Index (20°C)	1.405-1.409
Boiling Point	149-150°C
Acidity (CH ₃ COOH)	max. 0.02%
Water (KF)	max. 0.2%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No: 1110 , IMCO Class No.: 3 , Packing Group:III
HSN Code	
100 ml	29141990 (GST 18%)
500 ml	29141990 (GST 18%)

Available Packages

100 ml

500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com